

FOR OPTICIANS, OPTOMETRISTS AND CLINICS

CIOM

CIOM trial lens sets with 90, 158 or 232 lenses and the CIOM universal trial frame 120/T.

10 products · Prices on request: ask us for a quote

WEB SHOP

optinovatetech.com

WHATSAPP

+31 6 41861748

EMAIL

info@optinovatetech.com

Refraction

Trial sets

Trial lens sets for subjective refraction from CIOM: 90, 158 or 232 lenses in portable cases, beechwood cases or wooden trays, plus a universal trial frame and a concave/convex lens set.



CIOM

Trial Set 90 Lenses Portable case 100D

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- DIMENSIONS (W x D x H): 350x200x80 mm
- WEIGHT: Kg. 2
- This trial lens set offers a comprehensive selection of 90 lenses meticulously, the lightweight case...

Price on request: ask for a quote

[View online →](#)



CIOM

Trial Set Complete 232 Lenses, Beechwood Case 105/N

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- DIMENSIONS (W x D x H): 475x360x60 mm
- WEIGHT: Kg. 5
- Step into the realm of precision and versatility with our deluxe trial lens set, presenting a lavish ensemble...

Price on request: ask for a quote

[View online →](#)



CIOM

Trial Set Complete 232 Lenses, Lenses Only 108/L

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- The lenses, with their standard external diameter of 38 mm, can be easily inserted into any case or...

Price on request: ask for a quote

[View online →](#)



CIOM

Trial Set Complete 232 Lenses, Wooden Tray

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- DIMENSIONS (W x D x H): 445x310x45 mm
- WEIGHT: Kg. 4,7
- Featuring a complete array of 232 meticulously selected lenses, the wooden tray is thoughtfully designed to...

Price on request: ask for a quote

[View online →](#)



CIOM

Trial Set Complete 232 Lenses Portable Case 106/H

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- DIMENSIONS (W x D x H): 470x365x80 mm
- WEIGHT: Kg. 4,7
- Whether conducting routine eye examinations, fine-tuning contact lens fittings, or addressing complex visual...

Price on request: ask for a quote

[View online →](#)



CIOM

Trial Set Reduced to 158 Lenses, Lenses Only 103/L

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- The lenses, with their standard external diameter of 38 mm, can be easily inserted into any case or...

Price on request: ask for a quote

[View online →](#)



CIOM

Trial Set Reduced to 158 Lenses, Wooden Tray 103T

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- DIMENSIONS (W x D x H): 335x290x45 mm
- WEIGHT: Kg. 2
- Optimize your ophthalmic unit's efficiency with this compact trial lens set, featuring a reduced collection...

Price on request: ask for a quote

[View online →](#)



CIOM

Trial Set Reduced to 158 Lenses 101M

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- DIMENSIONS (W x D x H): 365x340x60 mm
- WEIGHT: Kg. 3,7
- This trial lens set offers a comprehensive selection of 158 lenses meticulously curated and presented in a...

Price on request: ask for a quote

[View online →](#)



CIOM

Trial Set Reduced to 158 Lenses Portable case 102P

- LENS: White optical glass, 22 mm diameter sphere / 22x18 mm cylinder
- RIM: Plastic, 38 mm diameter
- DIMENSIONS (W x D x H): 360x335x80 mm
- WEIGHT: Kg. 3,4
- Experience the perfect balance of portability and versatility with this medium-sized portable trial lens set...

Price on request: ask for a quote

[View online →](#)



CIOM

Universal Trial Frame 120/T

- Adjustable Pupillary Distance (PD): Ranges from 50 to 80 mm, accommodating various patient needs.
- Lens Holders: Equipped with five lens holders (three in front, two behind), allowing for comprehensive vision...
- 360° Lens Rotation: Ensures precise positioning of trial lenses.
- Adjustable Components: Features flexible temple ends and an adjustable nose rest for customized fitting.
- Lightweight Design: Weighing only 58 grams, it ensures patient comfort without compromising functionality.

Price on request: ask for a quote

[View online →](#)

How to order

1

Order online

Products with a price can be ordered and paid directly in our web shop, for example with iDEAL | Wero or a card.

2

Or request a quote

For several products, special wishes or products without a price: put them on your quote list and send it. You get a personal quote within 24 hours.

3

Questions?

WhatsApp +31 6 41861748 or email info@optinovatetech.com. We reply within 24 hours.

Prices and VAT

Prices on request: ask us for a quote. Quotes are in euro, excl. VAT. Businesses in another EU country with a valid VAT number pay 0% VAT when the goods are delivered to that country; outside the EU no Dutch VAT is charged.

Delivery

Delivery costs depend on the size of the product and the country. You see them at checkout, and all rates are on optinovatetech.com/delivery-and-vat/.

This catalogue was made automatically from our web shop on 20 September 2026. Products and prices can change; the web shop always shows the latest.

Optinovate Technologies · Turandotpad 7, 3194 WL Hoogvliet Rotterdam, the Netherlands

KVK 93248245 · VAT NL005010539B37 · info@optinovatetech.com · WhatsApp +31 6 41861748 · optinovatetech.com