

FOR OPTICIANS, OPTOMETRISTS AND CLINICS

Fundus examination

Fundus examination: aspheric lenses 20D, 78D and 90D, a Volk 3-mirror lens and the rechargeable direct ophthalmoscope TW-2411.

5 products · Prices on request: ask us for a quote

WEB SHOP

optinovatetech.com

WHATSAPP

+31 6 41861748

EMAIL

info@optinovatetech.com

Contents

Diagnostics

Ophthalmoscopes	1 product
Fundus examination	4 products

Optinovate Technologies is a Dutch company in Rotterdam that supplies eye-examination equipment to opticians, optometrists and clinics.

Every product name in this catalogue links to its page in our web shop, where you can order directly or put it on your quote list. We reply within 24 hours.

Diagnostics

Ophthalmoscopes

Rechargeable direct and indirect ophthalmoscopes for examining the fundus and retina: the direct ophthalmoscope TW-2411 and the indirect ophthalmoscopes TW-25C and TW-26.



Rechargeable Ophthalmoscope (Direct) TW-2411

- Specifications:
- 1. Illumination Source: 3.5V/2.8W halogen bulb
- 2. Recharging Unit Included
- 3. Input Voltage: 220V $\pm$ 10%
- 4. Input Power: 8VA

Price on request: ask for a quote

[View online →](#)

Fundus examination

Fundus examination: aspheric lenses 20D, 78D and 90D, a Volk 3-mirror lens and the rechargeable direct ophthalmoscope TW-2411.



Aspheric Lens 20D

- Image Magnification(X) 2.97
- Field of View 50 degree
- Working Distance(mm) 47
- Clear Aperture(mm) 48

Price on request: ask for a quote

[View online →](#)



Aspheric lens 78D

- Image Magnification(X) 0.95
- Field of View 88 degree
- Working Distance(mm) 7.0
- Clear Aperture(mm) 26

Price on request: ask for a quote

[View online →](#)



Aspheric lens 90D

- Image Magnification(X) 0.75
- Field of View 94 degree
- Working Distance(mm) 5.0
- Clear Aperture(mm) 18

Price on request: ask for a quote

[View online →](#)



Volk 3-Mirror Lens

Discover exceptional clarity and versatility with the Volk 3 Mirror Lens, an essential tool for comprehensive ophthalmic examinations. This innovative lens features three mirrors angled for precise visualization of the retina and peripheral structures. Crafted with premium materials, the Volk 3 Mirror Lens delivers high-resolution images and exceptional optical performance. Compatible with a wide...

Price on request: ask for a quote

[View online →](#)

How to order

1

Order online

Products with a price can be ordered and paid directly in our web shop, for example with iDEAL | Wero or a card.

2

Or request a quote

For several products, special wishes or products without a price: put them on your quote list and send it. You get a personal quote within 24 hours.

3

Questions?

WhatsApp +31 6 41861748 or email info@optinovatetech.com. We reply within 24 hours.

Prices and VAT

Prices on request: ask us for a quote. Quotes are in euro, excl. VAT. Businesses in another EU country with a valid VAT number pay 0% VAT when the goods are delivered to that country; outside the EU no Dutch VAT is charged.

Delivery

Delivery costs depend on the size of the product and the country. You see them at checkout, and all rates are on optinovatetech.com/delivery-and-vat/.

This catalogue was made automatically from our web shop on 20 September 2026. Products and prices can change; the web shop always shows the latest.

Optinovate Technologies · Turandotpad 7, 3194 WL Hoogvliet Rotterdam, the Netherlands

KVK 93248245 · VAT NL005010539B37 · info@optinovatetech.com · WhatsApp +31 6 41861748 · optinovatetech.com