



**OPTINOVATE**  
TECHNOLOGIES

CATALOGUE · SEPTEMBER 2026

FOR OPTICIANS, OPTOMETRISTS AND CLINICS

# Labomed

Labomed slit lamps (SL 20i, SLx 45, eVO300) and the Labomed Stella OPH microscope, supplied from Rotterdam.

4 products · Prices on request: ask us for a quote

WEB SHOP

**[optinovatetech.com](https://optinovatetech.com)**

WHATSAPP

**+31 6 41861748**

EMAIL

**[info@optinovatetech.com](mailto:info@optinovatetech.com)**

# Contents

## Slit lamps and microscopes

|                       |            |
|-----------------------|------------|
| Slit lamps            | 3 products |
| Operating microscopes | 1 product  |

**Optinovate Technologies** is a Dutch company in Rotterdam that supplies eye-examination equipment to opticians, optometrists and clinics.

Every product name in this catalogue links to its page in our web shop, where you can order directly or put it on your quote list. We reply within 24 hours.

# Slit lamps and microscopes

## Slit lamps

Buy Labomed LED slit lamps (slit lamp microscopes): the 2-step SL 20i for everyday screening, the 3-step SLx 45 up to 25× and the 5-step eVO300 up to 40× with a built-in photo port.



LABOMED

### eVO300 Slit Lamp, 5-Step, LED

- Equipped with an intelligent LED illumination system for enhanced visibility and clarity.
- Light Source: 5W LED, providing bright and consistent illumination.
- Binocular Tube: Features a 90° head for comfortable and precise viewing.
- Eyepiece: Wide Field (WF) 12.5X/18 mm with diopter adjustment range of  $\pm 5$ mm to accommodate different visual...
- Magnification Levels: Offers 5-step magnification options: 6.5X, 10X, 16X, 25X, and 40X, allowing for...

Price on request: ask for a quote

[View online →](#)



LABOMED

### SL 20i Slit Lamp, 2-Step, LED

- Galilean microscope with 2-step magnification: 10X and 16X
- Wide field eyepieces WF 12.5X/18mm with +/- 5 diopter adjustment
- Interpupillary distance range: 49 - 78mm
- Advanced LED illumination system (5W)
- Built-in filters: Red-free, cobalt blue, heat-absorbing

Price on request: ask for a quote

[View online →](#)



LABOMED

### SLx 45 Slit Lamp, 3-Step, LED

- Microscope Type: SLx 45 Slit Lamp with intelligent LED illumination system
- Binocular Tube: 90° head
- Eyepiece: WF 12.5X/18 mm, diopter adjustment  $\pm 5$ mm
- Magnification: 3-step options at 10X, 16X, and 25X
- Illumination: 5W LED

**Price on request: ask for a quote**

**[View online →](#)**

## Operating microscopes

Ophthalmic operating microscope: the Labomed Stella OPH microscope on a floor stand.



LABOMED

### Stella OPH Microscope, Floor Mount

- Model: Stella OPH microscope with intelligent LED illumination system
- Stand: Robust, vibration-free stand on cross base with four lockable castors
- Observation Tubes: Standard: 0-210° Binocular tiltable head, F=140mm
- Configured with Double beam splitter and Double Iris attachment for depth of field management
- Stereo Base: 24mm

**Price on request: ask for a quote**

**[View online →](#)**

# How to order

**1**

## Order online

Products with a price can be ordered and paid directly in our web shop, for example with iDEAL | Wero or a card.

**2**

## Or request a quote

For several products, special wishes or products without a price: put them on your quote list and send it. You get a personal quote within 24 hours.

**3**

## Questions?

WhatsApp +31 6 41861748 or email [info@optinovatetech.com](mailto:info@optinovatetech.com). We reply within 24 hours.

## Prices and VAT

Prices on request: ask us for a quote. Quotes are in euro, excl. VAT. Businesses in another EU country with a valid VAT number pay 0% VAT when the goods are delivered to that country; outside the EU no Dutch VAT is charged.

## Delivery

Delivery costs depend on the size of the product and the country. You see them at checkout, and all rates are on [optinovatetech.com/delivery-and-vat/](https://optinovatetech.com/delivery-and-vat/).

This catalogue was made automatically from our web shop on 20 September 2026. Products and prices can change; the web shop always shows the latest.

---

**Optinovate Technologies** · Turandotpad 7, 3194 WL Hoogvliet Rotterdam, the Netherlands

KVK 93248245 · VAT NL005010539B37 · [info@optinovatetech.com](mailto:info@optinovatetech.com) · WhatsApp +31 6 41861748 · [optinovatetech.com](https://optinovatetech.com)