



OPTINOVATE
TECHNOLOGIES

CATALOGUE · SEPTEMBER 2026

FOR OPTICIANS, OPTOMETRISTS AND CLINICS

Tonometers

Tonometers for measuring intraocular pressure (IOP): the applanation tonometer TW-2417R/T, the indentation tonometer TW-2470 and the rebound tonometer IN V033.

3 products · Prices on request: ask us for a quote

WEB SHOP

optinovatetech.com

WHATSAPP

+31 6 41861748

EMAIL

info@optinovatetech.com

Diagnostics

Tonometers

Tonometers for measuring intraocular pressure (IOP): the applanation tonometer TW-2417R/T, the indentation tonometer TW-2470 and the rebound tonometer IN V033.



Applanation Tonometer TW-2417R/T

- Features-
- The product is designed on the principle basis of Goldman tonometer.
- Three types are available: R and T type are suitable for Haag Streit type slit lamp, While X type is suitable...

Price on request: ask for a quote

[View online →](#)



Indentation Tonometer TW-2470

Tonometer is the mostly used device in hospital at present. Each division on the scale corresponds to an indentation of 1/20 mm on the cornea.

Price on request: ask for a quote

[View online →](#)



Rebound Tonometer IN V033

- Rebound Tonometer IN-V033
- Accurate. Portable. Comfortable.
- Key Features:
- Measuring Range: -3 mmHg to 70 mmHg
- Precision: ± 1.5 mmHg for intraocular pressures ≤ 25 mmHg

Price on request: ask for a quote

[View online →](#)

How to order

1

Order online

Products with a price can be ordered and paid directly in our web shop, for example with iDEAL | Wero or a card.

2

Or request a quote

For several products, special wishes or products without a price: put them on your quote list and send it. You get a personal quote within 24 hours.

3

Questions?

WhatsApp +31 6 41861748 or email info@optinovatetech.com. We reply within 24 hours.

Prices and VAT

Prices on request: ask us for a quote. Quotes are in euro, excl. VAT. Businesses in another EU country with a valid VAT number pay 0% VAT when the goods are delivered to that country; outside the EU no Dutch VAT is charged.

Delivery

Delivery costs depend on the size of the product and the country. You see them at checkout, and all rates are on optinovatetech.com/delivery-and-vat/.

This catalogue was made automatically from our web shop on 20 September 2026. Products and prices can change; the web shop always shows the latest.

Optinovate Technologies · Turandotpad 7, 3194 WL Hoogvliet Rotterdam, the Netherlands

KVK 93248245 · VAT NL005010539B37 · info@optinovatetech.com · WhatsApp +31 6 41861748 · optinovatetech.com